



Vision Tools Workspace

Integrating the Healthcare Enterprise (IHE) Integration Statement

Millensys Contacts

Main Office

Address: 36 Ahmed El Zommor St. Nasr City – Cairo

Tel / Fax: +20 2 22713860 - +20 2 22754867 - +20 2 22879403

Web site

www.MillenSys.com

Email

Sales : sales@millensys.com

Support : support@millensys.com

All material contained in this document is copyright Millensys.

Millensys and Vision Tools Logo are trademarks of Millensys.

Microsoft and Windows are trademarks of Microsoft Corporation in the United States and other countries.

The information presented in this document is correct at time of going to press. However, it is subject to change without notice.

This document may not be photocopied, reproduced or translated into another language, in whole or in part, without the consent of Millensys.

MillenSys Workspace IHE. 

Document Information

| | | |
|---|------------------------|--|
| Author | Eng.Ahmed Hussein Riad | Tel: +20 2 29287061 E-Mail: a_hussein@millensys.com |
| Editing & layout | Eng.Hossam Rady | Tel: +20 2 29287061 E-Mail: h_rady@millensys.com |
| For questions regarding the technical content of this document, please contact: | Eng.Ahmed Hussein Riad | Tel: +20 2 29287061 E-Mail: a_hussein@millensys.com |

Table 1: Document Information Table

| IHE Integration Statement | | Date Published 01 February, 2017 |
|---|----------------------------------|--|
| Vendor | Product Name | Product Version |
| Millensys Systems | Vision Tools Workspace | 5.X |
| This product implements all transaction required in the IHE framework to support the IHE Integration Profile, Actors and Options listed below | | |
| Integration Profile Implements | Actors Implemented | Options Implemented |
| Consistent Time (CT) | Time Client | |
| Scheduled Workflow (SWF) | Image Manager\Image Archive | Assisted Acquisition Protocol Setting |
| | Order Filler | |
| | Performed Procedure Step Manager | |
| Access To Radiology Information System (ARI) | Image Manager\Image Archive | |
| Consistent Presentation of Images (CPI) | Image Manager\Image Archive | |

| Integration Profile Implements | Actors Implemented | Options Implemented |
|--|---|---------------------|
| Mammography Image (MAMMO) | Image Manager\Image Archive | |
| | Image Display | |
| Digital Breast Tomosynthesis (DBT) | Image Manager\Image Archive | |
| NM Image (NMI) | Image Manager\Image Archive | |
| Echocardiography | Image Manager\Image Archive | |
| Patient Information Reconciliation (PIR) | Image Manager\Image Archive | |
| | Performed Procedure Step Manager | |
| Evidence Documents (ED) | Image Manager\Image Archive | |
| Basic Image Review (BIR) | Image Manager\Image Archive | |
| Key Images Notes (KIN) | Image Manager\Image Archive/Image Archive | |
| | Evidence Creator | |
| | Image Display | |
| Portable Data for Imaging (PDI) | Portable Media Importer | |
| Report Workflow (RWF) | Report Creator | |
| | Report Reader | |
| | Report Manager | |
| | Department System Scheduler\Order Filler | |

| Integration Profile Implements | Actors Implemented | Options Implemented |
|--|--|---------------------|
| Audit Trail and node Authentication (ATNA) | Image Manager\Image Archive | |
| | Image Display | |
| | Department System Scheduler\Order Filler | |
| Presentation of Grouped Procedures (PGP) | Department System Scheduler\Order Filler | |
| | Image Manager\Image Archive | |
| | Performed Procedure Step Manager | |

References

<http://www.rsna.org/IHE>

Vendor's IHE information

http://www.millensys.com/docs/workspace_ihe.pdf

Vendor's DICOM information

http://www.millensys.com/docs/workspace_dicom.pdf

General Information on IHE

North America : <http://www.rsna.org/IHE>

Europe : <http://www.ihe-europe.org>

Japan : <http://www.jira-net.or.jp/ihe-j>